

BACK EXTENSION

Tools for the improvement of active rehabilitation




BACK EXTENSION

- TOTAL WEIGHT 150 kg
- LENGTH 125 cm
- WIDTH 100 cm
- HEIGHT 206 cm
- WEIGHT STACK 50 kg
- optional 75 kg

Complies to CE-MDD regulations Class 1 (with measurement function).

The fully adjustable inclination of the pelvis is one of the unique features of the Monitored Rehab Systems Back Extension. In combination with the mechanical ROM limiter and the MRS Software the Back Extension is ideally suited for a wide range of therapeutic goals. From stabilizing the lower back to mobilizing over a wider ROM the Back Extension ensures maximal adaptation to your therapeutic intervention. Visual feedback stimulates patient effort and helps the clinician to track progress.



ISOMETRIC GATE

A dynamic test with an isometric 'hold', this exercise guides patients as they maneuver the red ball between the gates in the bars.



CONTROLLED ROUTE

Concentric and eccentric movements shift the circle up and down, activating muscle control and invigorating coordinate abilities.



MOTOR CONTROLE

Eccentric and Concentric exercise settings allows clinician to control speeds both concentrically and eccentrically affecting motor control and joint stability.



MATHEMATICS DUAL TASKS

The Dual Tasking option incorporates math problems and movement of various muscle groups to complete a task.



SKI GAME

Full 'Super G' downhill. A very challenging exercise.

To schedule a demonstration, contact your local dealer.

mrsystems.nl